

PATHFAST cTnI Meets the Criteria for High-Sensitivity Troponin Assays

E. Spanuth¹, R. Thomae², E. Giannitsis³

¹DIAneering®, Diagnostics Engineering & Research GmbH, Heidelberg, Germany, ²Mitsubishi Chemical Europe, Düsseldorf, Germany, ³Dept. of Medicine III, Div. of Cardiology, University of Heidelberg, Heidelberg, Germany

Background

Contemporary available point-of-care assays for troponin cannot provide the analytical and diagnostic accuracy required for “highly-sensitivity” assays by international guidelines. The PATHFAST cTnI POC assay has shown promising analytical validity.

Aim of the study

To examine whether the PATHFASTcTnI could be classified “highly-sensitive”.

Methods

Control group: 119 healthy individuals (60 men and 59 women, 21-69 years old, median 42 years).

Patient group: cTnI (PATHFAST) and cTnT (cobas® hs-cTnT) were compared in 181 patients admitted to the chest pain unit at presentation, 3 and 6 hours later. The results were related to the discharge diagnoses.

Results

Tab. 1: PATHFAST cTnI (ng/L) in healthy individuals (control group)

	n	Mean	IQR	Lowest value	Highest value
All	119	2.06	0.87-1.86	0.37	17.2
Men	60	2.82	1.15-2.64	0.53	17.2
Women	59	1.12	0.70-1.25	0.37	11.2

p < 0.0001

> 99th percentile:

16 ng/L (CLSI C28-A3)
manufacturer recommended: 20.0 ng/L

> Quantification:

66.4% (> LoD (1.0 ng/L) and < 99th percentile):
(in 79 of 119 healthy individuals of the control group)

> Imprecision:

cTnI (ng/L)	CV (NCCLS)
2.0	20%
3.0	10%
20.0	5%

> Imprecision Profile:




Fig. 1: Imprecision Profile according to NCCLS

A: at 1.93, 3.13, 8.29, 15.32 and 19.66 ng/L

B: CVs from 3-fold determination of controls, n=119

Tab. 2: PATHFAST cTnI in patients with chest pain at presentation (0 h), and after 3 (3 h) and 6 hours (6 h)

	NSTEMI, n=72 Median (IQR)	UAP and NCCP* Median (IQR)
0 h	46 (34-76)	2.0 (1.0-4.0)
3 h	166 (52-759)	2.0 (2.0-5.0)
6 h	399 (85-1100)	3.0 (2.0-6.0)

*UAP (unstable angina pectoris), n=80; NCCP (non cardiac chest pain), n=29

Tab. 3: Results of ROC analysis for diagnosis of NSTEMI

	PATHFAST cTnI			hs-cTnT (cobas)		
	AUC	SENS (%)	SPEC (%)	AUC	SENS (%)	SPEC (%)
0 h	0.919	83	93	0.923	77	95
3 h	0.962	93	96	0.964	92	97
6 h	0.958	91	94	0.969	94	97

Absolute change of cTnI:

	cTnI < URL at admission (n=16)			cTnI ≥ URL at admission (n=56)		
	Median	Change		Median	Change	
0h	0,0050	(pg/ml)	(%)	0.099	(pg/ml)	(%)
3h	0,0405	0.036	710	0.328	0.229	231
6h	0,141	0.101	248	0.524	0.196	59.8




Fig. 2: ROC curve of absolute change for NSTEMI

A: between baseline and 3 hours; B: between baseline and 6 hours

Conclusion

PATHFAST cTnI meets the criteria for high-sensitive cTn assays:

- > Imprecision CV at the 99th percentile: 5% (<10%)
- > Quantification in healthy individuals: 66.4% (>50%)
- > Complete concordance with cobas hs-cTnT for the diagnosis of NSTEMI
- > Absolute changes of cTnI values over time in NSTEMI patients > 30%