



Antibodies against Phospholipids

for Diagnosis of the Anti-Phospholipid-Syndrome

Thomas Büttner

Product Manager



- Situated near Berlin
- Founded 1992 Medipan, 2002 GA
- About 80 employees
- Manufacturer of test systems for autoimmune diagnostics

Anti-Phospholipid-Syndrome

Definition

The **Anti-Phospholipid Syndrome** (APS) is one of the most common autoimmune diseases.

The disease is characterized by the appearance of specific **antibodies against various phospholipids** (e.g. cardiolipins) or **phospholipid-binding proteins** (e.g. prothrombin or β 2-glycoprotein I).

These antibodies cause an **increased tendency to coagulate (hypercoagulability)** and lead to increased **thrombosis**.



Anti-Phospholipid-Syndrome

Definition

Primary APS (PAPS) Independent from other diseases

Secodary APS (SAPS) occurs in the context of other, mainly autoimmune, diseases (in particular SLE, RA, SjS) or tumours

Antibodies against Phospholipids for the Diagnosis of the Anti-Phospholipid-Syndrome



Anti-Phospholipid-Syndrome

Epidemiologie

Worldwide, 2 - 5% of the population are affected, primarily women.

Anti-Phospholipid-Syndrome

Clinical Symptoms

The clinical picture is characterized by

- an **increased occurrence of thromboses and embolisms** as well as pulmonary embolisms, heart attacks, strokes, renal infarctions, etc.,
- **bleeding and ulcers of the skin** (purpura),
- **miscarriages and early abortions**
- **most severe impact:** catastrophic APS (CAPS), extensive thromboses in small vessels that supply multiple organs


Anti-Phospholipid-Syndrome

Therapy

The treatment of asymptomatic patients consists of **thrombosis prophylaxis**, e.g. with acetylsalicylic acid or hydroxychloroquine.

After a thrombosis, the more effective **anti-coagulant therapy** is the administration of vitamin K antagonists such as phenprocoumon (Marcumar®, Falithrom®).


If necessary, the lifelong therapy is supported by the administration of **immunosuppressants**.



Anti-Phospholipid-Syndrome

Classification of Disease

Sapporo criteria, 1999



Anti-Phospholipid-Syndrome

Classification of Disease

Sapporo criteria, 1999

- Presence of clinical symptoms
 - One or more episodes of arterial, venous or small-vessel thrombosis in any tissue or organ; confirmation by imaging or histopathological procedures
 - Pregnancy diseases (stillbirths, premature births, spontaneous abortions)
- Laboratory parameter
 - Cardiolipin antibodies of IgG and/or IgM type, medium or high titre measured by a standardised, β 2-glycoprotein-dependent ELISA
 - Positive Lupus anticoagulant in plasma,
- detection at two or more points in time at least 6 weeks apart;

Anti-Phospholipid-Syndrome

Classification of Disease

Sapporo criteria, 1999

Revised criteria Sydney, 2004

Journal of Thrombosis and Haemostasis, 4: 295–306

SPECIAL ARTICLE

International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS)

S. MIYAKIS,* M. D. LOCKSHIN,† T. ATSUMI,‡ D. W. BRANCH,§ R. L. BREY,¶ R. CERVERA,**
R. H. W. M. DERKSEN,†† P. G. DE GROOT,†† T. KOIKE,‡ P. L. MERONI,‡‡ G. REBER,§§
Y. SHOENFELD,¶¶ A. TINCANI,*** P. G. VLACHOYIANNOPOULOS††† and S. A. KRILIS*

*St George Hospital, University of New South Wales, Sydney, Australia; †Hospital for Special Surgery, Cornell Medical Center, New York, NY, USA; ‡Hokkaido University, Sapporo, Japan; §University of Utah Health Sciences Center, Salt Lake City, UT; ¶University of Texas Health Science Center, San Antonio, TX, USA; **Hospital Clinic, Barcelona, Spain; ††University Medical Center, Utrecht, The Netherlands; ‡‡Istituto Auxologico Italiano, University of Milan, Milan, Italy; §§University Hospital, Geneva, Switzerland; ¶¶Sheba Medical Center, Tel-Hashomer and Tel Aviv University, Israel; ***Spedali Civili, University of Brescia, Italy; and †††Department of Pathophysiology, University of Athens, Greece

To cite this article: Miyakis S, Lockshin MD, Atsumi T, Branch DW, Brey RL, Cervera R, Derksen RHWM, de Groot PG, Koike T, Meroni PL, Reber G, Shoenfeld Y, Tincani A, Vlachoyiannopoulos PG, Krilis SA. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost* 2006; 4: 295–306.

Summary. New clinical, laboratory and experimental insights, since the 1999 publication of the Sapporo preliminary classification criteria for antiphospholipid syndrome (APS), had been addressed at a workshop in Sydney, Australia, before the

for APS. Members of the workshop panel included all of the authors and the individuals listed in the Appendix.

Some of the authors presented the current evidence in their area of expertise (see Addendum) providing relevant literature on predictors of outcome, risk factors, associations between

Anti-Phospholipid-Syndrome

Classification of Disease

Sapporo criteria, 1999

Revised criteria Sydney, 2004

- Additional Laboratory parameter
 - Positive detection of antibodies to β 2GPI of the IgG and/or IgM type at least twice in serum or plasma with a titre level above the 99th percentile is to be classified as a laboratory analytical criterion for APS using standardised ELISA.
- Check-up after initial antibody detection should be carried out at least every 12 weeks (instead of the previous 6 weeks)

Anti-Phospholipid-Syndrome

Classification of Disease



Detection of antibodies necessary for classification of Anti-Phospholipid syndrome beside clinical symptoms

Differentiation from other thrombophilic events !

Anti-Phospholipid-Syndrome

History of APS Diagnostics

1906: reaction of sera from Syphilis patients with preparation from bovine heart (Cardiolipin) – complement fixation



Anti-Phospholipid-Syndrome

History of APS Diagnostics

1946: Further development of the Wassermann test by Harris, Rosenberg and Riedel
Venereal disease Research Laboratory (VDRL Test)

➡ formation of anti-Cardiolipin Antibodies following the Treponema infection

Later: false positive VDRL tests in patients with

- Collagenoses (especially Systemic Lupus erythematosus)
- other infections: Lepra, Mononucleosis infectiosa, Virus hepatitis

Anti-Phospholipid-Syndrome

History of APS Diagnostics

1980s: solid-phase assays for the detection of anti-phospholipid antibodies

Volume 322, Issue 8361, 26 November 1983, Pages 1211-1214

ANTICARDIOLIPIN ANTIBODIES: DETECTION BY RADIOIMMUNOASSAY AND ASSOCIATION WITH THROMBOSIS IN SYSTEMIC LUPUS ERYTHEMATOSUS

E.N. Harris, M.L. Boey, C.G. Mackworth-Young, A.E. Gharavi, B.M. Patel, S. Loizou, G.R.V. Hughes

Show more ▾

+ Add to Mendeley  Share  Cite

[https://doi.org/10.1016/50140-6736\(83\)91267-9](https://doi.org/10.1016/50140-6736(83)91267-9)

[Get rights and content](#)

Abstract

A new solid-phase radioimmunoassay for the detection of anticardiolipin antibodies is 200-400 times more sensitive than the precipitation method used in the Venereal Disease Reference Laboratory test. 61% of serum samples from patients with systemic lupus erythematosus (SLE) had high levels of anticardiolipin antibodies of at least one immunoglobulin class. There were strong correlations between raised anticardiolipin

Clin. exp. Immunol. (1984) **56**, 193-199.

Anti-phospholipid antibodies and biological false positive serological test for syphilis in patients with systemic lupus erythematosus

T. KOIKE, M. SUEISHI, H. FUNAKI, H. TOMIOKA & S. YOSHIDA *Second Department of Internal Medicine, School of Medicine, Chiba University, Chiba, Japan*

(Accepted for publication 31 October 1983)

SUMMARY

Utilizing a newly developed solid phase enzyme immunoassay for the detection of anti-phospholipid antibodies, we found that 41.7% of sera from patients with active systemic lupus erythematosus (SLE) had the anti-cardiolipin antibody and 22.7% of sera were positive for the anti-phosphatidylinositol (PI) antibody. But antibodies to lecithin, cephalin and sphingomyelin were rarely found in this assay. We also examined antibodies

Anti-Phospholipid-Syndrome

History of APS Diagnostics

1990s: cofactor for anti-phospholipid antibodies in autoimmune diseases

Proc. Natl. Acad. Sci. USA
Vol. 87, pp. 4120–4124, June 1990
Medical Sciences

Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: β_2 -Glycoprotein I (apolipoprotein H)

(anti-cardiolipin antibodies/lipid-protein complex/thrombosis)

H. PATRICK McNEIL*, RICHARD J. SIMPSON†, COLIN N. CHESTERMAN*, AND STEVEN A. KRILIS*‡

*University of New South Wales, School of Medicine, Saint George Hospital, Kogarah, 2217 Australia; and †Joint Protein Structure Laboratories, The Ludwig Institute for Cancer Research and The Walter and Eliza Hall Institute for Medical Research, Parkville, 3050 Australia

Communicated by K. Frank Austen, March 15, 1990

ABSTRACT Anti-phospholipid (aPL) antibodies that exhibit binding in cardiolipin (CL) ELISA can be purified to >95% purity by sequential phospholipid affinity and ion-exchange chromatography. However, these highly purified aPL antibodies do not bind to the CL antigen when assayed by a modified CL ELISA in which the blocking agent does not contain bovine serum, nor do they bind to phospholipid affinity

We have previously noted that when aCL antibody-containing fractions derived from ion-exchange chromatography of plasma were applied to phosphatidylserine or CL affinity columns, there was no binding of the antibody despite the fact that when plasma containing these antibodies was applied to these columns, aCL antibodies could be purified. This suggested that there was a cofactor also present in

> *J Immunol.* 1992 Jun 15;148(12):3885-91.

Heterogeneity of anticardiolipin antibodies defined by the anticardiolipin cofactor

E Matsuura¹, Y Igarashi, M Fujimoto, K Ichikawa, T Suzuki, T Sumida, T Yasuda, T Koike

Affiliations + expand

PMID: 1602135

Abstract

Anticardiolipin antibodies (aCL) found in sera from patients with SLE react with cardiolipin (CL) in the presence of a 50-kDa serum cofactor. The cofactor, which was identified to be beta 2-glycoprotein I by sequencing the N-terminal amino acids, not only enhances CL binding by antibodies in SLE but also depresses it by antibodies associated with syphilis. Cofactor-dependent binding of aCL in SLE to solid phase CL was competitively inhibited by the simultaneous addition of fluid phase CL but was unaffected by either prior or simultaneous addition of a high excess of the cofactor. Binding of aCL in

Anti-Phospholipid-Syndrome

History of APS Diagnostics

1990s: cofactor for anti-phospholipid antibodies in autoimmune diseases

> [Ann Ital Med Int.](#) 1993 Jul-Sep;8(3):171-4.

Anti-cardiolipin antibodies in HIV infection are true antiphospholipids not associated with antiphospholipid syndrome

M Falco ¹, A Sorrenti, R Priori, F L Luan, V Pittoni, M G Agresti, G Valesini

Affiliations + expand

PMID: 8217481

Abstract

The aim of the present study was to evaluate the fine specificity of anticardiolipin (aCL) antibodies detectable in the sera of patients with HIV infection. aCL are generally associated with thrombotic events in autoimmune diseases. A solid phase ELISA which discriminates between aCL binding to phospholipids and aCL binding to phospholipid/beta 2-glycoprotein I (cofactor) complex was employed. Thirty-nine HIV and 20 aCL positive systemic lupus erythematosus (SLE) sera were examined. In HIV sera, reduced binding to phospholipid was seen if cofactor was added. On the contrary, in SLE-sera the cofactor improved aCL binding. No thrombotic events were recorded in HIV

Differences between antibodies between autoimmunity and infection !

Anti-Phospholipid-Syndrome

History of APS Diagnostics

1990s: cofactor for anti-phospholipid antibodies in autoimmune diseases

[J Exp Med.](#) 1994 Feb 1;179(2):457-62. doi: 10.1084/jem.179.2.457.

Anticardiolipin antibodies recognize beta 2-glycoprotein I structure altered by interacting with an oxygen modified solid phase surface

E Matsuura ¹, Y Igarashi, T Yasuda, D A Triplett, T Koike

Affiliations + expand

PMID: 7507506 PMCID: PMC2191370 DOI: 10.1084/jem.179.2.457

[Free PMC article](#)

Abstract

Anticardiolipin antibodies (aCL) derived from the sera of individuals exhibiting the antiphospholipid syndrome (APS) directly bind to beta 2-glycoprotein I (beta 2-GPI), which is adsorbed to an oxidized polystyrene surface. Oxygen atoms were introduced on a polystyrene surface by irradiation with electron or gamma-ray radiation. X-ray photoelectron spectroscopy revealed the irradiated surfaces

No phospholipids necessary for antibody determination !



Anti-Phospholipid-Antibodies

3 Types of antibodies under the same name!

- Antibodies to phospholipids
- Antibodies to phospholipid.cofactor complexes
- Antibodies to cofactor proteins

Standardisation ?

Anti-Phospholipid-Syndrome

Standardisation ?

E.J. Favaloro (Dept. of Haematology, Institute of Clinical Pathology and Medical Research, Western Sydney Area Health Service, Westmead, NSW, Australia)

RCPA Immunology Newsletters, 5, 2, 2003, pp. 3-6

Testing of 7 commercial Anti-Cardiolipin Assays in 56 laboratories (RCPA Quality Assurance Program)




What antibody are we talking about ?

Phospholipids

Structure

Phospholipids, e.g. phosphatidylcholine (lecithin), phosphatidylethanolamine, phosphatidylserine or phosphatidylinositol, are **polar, phosphorylated lipids** in the cell membrane.




Phosphatidylserin

Cardiolipins

Structure

Cardiolipins are phospholipids with four fatty acid residues.

They are the only phospholipids synthesized in the **mitochondria**.



Cardiolipin

Antibodies against Phospholipids

Cofactor β 2-glycoprotein I


Plasma protein: 200 mg/L

Molecular weight: about 50 kDa

may play a role in coagulation, bind to platelets, mitochondria and negatively charged substances such as heparin, DNA, dextran sulphate and negatively charged phospholipids

β 2-Glycoprotein I

326 amino acids organized in 5 domains (DI – DV)



opening of circular structure by binding to negatively charged phospholipid

hidden epitope on domain I available for antibodies

Ruben E, Planer W: J. Biol. Chem. (2020) 295(31) 10794–10806

Antibodies against Phospholipids

Standardisation?

UK NEQAS

Immunology, Immunochemistry & Allergy

Sample 234-1 was a single donor serum from a patient with an elevated concentration of IgM anti-phospholipid antibodies diluted with pooled normal human serum.

Sample	Analyte	Target Response
234-1	ACA IgG	Negative
	ACA IgM	Positive
	Anti-B2-Glycoprotein IgG	Equivocal
	Anti-B2-Glycoprotein IgM	Positive
	Anti-Phosphatidyl Serine IgG	Negative
	ACA (Screen Only)	Positive

Antibodies against Phospholipids

Standardisation?

Sample	Analyte	Target Response
234-1	ACA IgG	Negative
	ACA IgM	Positive
	Anti-B2-Glycoprotein IgG	Equivocal
	Anti-B2-Glycoprotein IgM	Positive
	Anti-Phosphatidyl Serine IgG	Negative
	ACA (Screen Only)	Positive

Anti Phospholipid Antibodies - Sample 234-1						
Distribution :	234	July 2023			Participant	90323
ACA IgG						
Your Results						
ACA IgG:	Negative					
Quantitation:	2.1					
Unit:	Units/ml					
Method:	Enzyme Immunoassay/ELISA					
Manufacturer:	Generic Assays					
	Returns	Positive	Negative	Mean	Range	
Other ...						
Other ...	4	1	3	2	0	- 4
Other ...	4	1	3	2	0	- 4
Units/ml						
Chemiluminescence	117	0	117	8	0	- 18
A. Menarini - Zenit RA	22	0	22	1	0	- 6
HOB-BioCLIA	2	0	2	6	5	- 7
IDS-ISYS	11	0	11	1	0	- 4
Inova QUANTA Flash	77	0	77	11	0	- 18
Other ...	5	0	5	10	6	- 15
Enzyme Immunoassay/ELISA	283	0	283	4	0	- 27
A. Menarini - Zenit EIA	2	0	2	5	0	- 10
AESKU	16	0	16	5	1	- 10
Biosystems	1	0	1	2	2	- 2
CLS AUTOZYME	3	0	3	8	6	- 10
Diesse	1	0	1	4	4	- 4
Euroimmun	9	0	9	3	2	- 5
Generic Assays	2	0	2	2	2	- 2
Human-Intec	2	0	2	17	6	- 27
Immunoconcepts	1	0	1	6	6	- 6
In House	1	0	1	4	4	- 4
Inova	21	0	21	4	1	- 9
ORGenTec	24	0	24	3	1	- 4
ORGenTec Alegria	8	0	8	4	4	- 5
Other ...	4	0	4	6	6	- 7
Phadia 100 EliA	2	0	2	3	3	- 3
Phadia 1000/2500 EliA	10	0	10	3	3	- 4
Phadia 250 EliA	176	0	176	3	2	- 12
Multiple Particle Based Assays	36	35	1	30	9	- 36
Bio-Rad Bioplex 2200	35	35	0	30	23	- 36
THERADIAG (BMD)	1	0	1	9	9	- 9

Antibodies against Phospholipids

Standardisation?

Sample	Analyte	Target Response
234-1	ACA IgG	Negative
	ACA IgM	Positive
	Anti-B2-Glycoprotein IgG	Equivocal
	Anti-B2-Glycoprotein IgM	Positive
	Anti-Phosphatidyl Serine IgG	Negative
	ACA (Screen Only)	Positive

Anti Phospholipid Antibodies - Sample 234-1

Distribution :	234	July 2023	Participant	90323
----------------	-----	-----------	-------------	-------

Anti-B2 Glycoprotein IgG

Your Results
Anti-B2 Glycoprotein IgG: Negative
Quantitation: 4.4
Unit: Units/ml
Method: Enzyme Immunoassay/ELISA
Manufacturer: Generic Assays

	Returns	Positive	Negative	Mean	Range		
Units/ml							
Chemiluminescence	112	74	38	25	0	-	60
A. Menarini - Zenit RA	22	0	22	1	0	-	3
HOB-BioCLIA	1	0	1	8	8	-	8
IDS-iSYS	10	0	10	1	0	-	1
Inova QUANTA Flash	74	70	4	34	0	-	60
Other ...	5	4	1	27	5	-	45
Enzyme Immunoassay/ELISA	246	1	245	4	0	-	26
A. Menarini - Zenit EIA	1	0	1	1	1	-	1
AESKU	15	0	15	4	0	-	8
Corgenix	1	0	1	4	4	-	4
Diesse	1	0	1	5	5	-	5
Euro Diagnostica	1	0	1	4	4	-	4
Euroimmun	6	0	6	3	2	-	3
Generic Assays	2	0	2	2	0	-	4
Human-Imtec	2	0	2	4	3	-	5
In House	1	0	1	1	1	-	1
Inova	16	0	16	3	0	-	5
ORGenTec	19	0	19	2	1	-	4
ORGenTec Alegria	6	0	6	5	4	-	6
Other ...	3	0	3	3	0	-	6
Phadia 100 EliA	2	0	2	3	3	-	4
Phadia 1000/2500 EliA	13	0	13	5	1	-	6
Phadia 250 EliA	157	1	156	4	1	-	26
Multiple Particle Based Assays	35	35	0	26	14	-	30
Bio-Rad Bioplex 2200	34	34	0	26	22	-	30
THERADIAG (BMD)	1	1	0	14	14	-	14

Antibodies against Phospholipids

Standardisation?


Cardiolipin Antibodies

accepted protocol: cofactor dependent antigen

β 2-glycoprotein I domain I specific assays?

QUANTA Flash® β 2GP1 IgG

QUANTA Flash® β 2GP1 Domain 1




P. de Groot, Utrecht

Anti-Phospholipid-Syndrome

Diagnosis

Diagnosis is based on **clinical symptoms** and **laboratory findings**.

In laboratory diagnostics, detection is primarily carried out by the determination of **antibodies against phospholipids** (e.g. against cardiolipin) and **associated proteins** (e.g. β 2-glycoprotein I).



Anti-Cardiolipin

Assays from Generic Assays (Elisa)

The **Anti-Cardiolipin** is a ELISA immunoassay for the quantitative determination of IgG and/or IgM antibodies against cardiolipin- β 2GPI in human serum.



Anti-Cardiolipin screen

Assays from Generic Assays (Elisa)

The **Anti-Cardiolipin screen** is a ELISA immunoassay for the semi-quantitative determination of IgG/IgM/IgA antibodies against cardiolipin- β 2GPI in human serum.



Product Portfolio

ELISA for the Diagnosis of Anti-Phospholipid-Syndrome

IgA + IgG + IgM Detection

Anti-Cardiolipin Screen

Anti- β 2-GP-I Screen

IgG or IgM Detection

Anti-Cardiolipin

Anti- β 2-GP-I


Anti-Phosphatidyl-Serin

Anti-Phospholipid Screen


Anti-Phospholipid Dot (IgG, IgM)

Multi-parameter Assay from GA Generic Assays

The **Anti-Phospholipid Dot** is a line immunoassay (Elisa principle) for the qualitative determination of IgG / IgM antibodies against 7 phospholipids and 3 serum proteins in human serum.




Anti-Phospholipid Dot (IgG / IgM)



Prof. René Louis Humbel, Luxemburg

Anti-Phospholipid Dot (IgG / IgM)

Anti-Phospholipid 10 Dot (REF 5012)
Patients with manifest APS:
IgG results:



- often reaction of all anionic phospholipids and β 2-GPI

DotDiver Anti-Phospholipid:

VDRL pos. sera (Treponema infection):



no reaction with β 2-GPI



Anti-Phospholipid Dot (IgG / IgM)

Lupus (2020) 0, 1–9
journals.sagepub.com/home/lup

PAPER

Profiling of non-criteria antiphospholipid antibodies in patients with SLE: differentiation of thrombotic SLE patients and risk of recurrence of thrombosis

O Tkachenko¹, S Lapin¹, A Mazing¹, V Emanuel¹, E Belolipetskaia², I Beliaeva², V Myachikova³, A Maslyansky³, P Schierack^{4,5} and D Roggenbuck^{4,5}

¹Center for Molecular Medicine, First Pavlov State Medical University of Saint Petersburg, Saint Petersburg, Russian Federation; ²North-Western State Medical University named after II Mechnikov, Saint Petersburg, Russian Federation; ³Rheumatology Department, VA Almazov North-West Federal Medical Research Center, Saint Petersburg, Russian Federation; ⁴Faculty Environment and Natural Sciences, Brandenburg University of Technology Cottbus-Senftenberg, Senftenberg, Germany; and ⁵Faculty of Health Sciences, Joint Faculty of the Brandenburg University of Technology Cottbus – Senftenberg, the Brandenburg Medical School Theodor Fontane and the University of Potsdam, Germany

To reveal the clinical significance of criteria and non-criteria antiphospholipid antibodies detected by line immunoassay in comparison with ELISA, systemic lupus erythematosus patients with and without thrombotic events were investigated. Thus, 107 systemic lupus erythematosus patients (48% with deep vein thrombosis or/and arterial thrombosis) and 120 healthy donors were enrolled. Serum antiphospholipid antibodies were detected by ELISA (Orgentec Diagnostika, Germany) and line immunoassay (GA Generic Assays, Germany). Lupus anticoagulant and IgG to cardiolipin and β 2GPI but not IgM as well as triple positivity by ELISA and line immunoassay were linked with thrombosis in systemic

Why checking for 10 parameters ?

Diagnostic advantages of detection of multiple parameters

Anti-Phospholipid Dot (IgG / IgM)

Why checking for 10 parameters ?

Diagnostic advantages of detection of multiple parameters




Figure 2 Comparison of quantitative anti cardiolipin (aCL) IgM (a) and anti-beta 2 glycoprotein I (aB2GPI) IgM (b), aCL IgG (c), aB2GPI IgG (d) analysis by line immunoassay (LIA) in systemic lupus erythematosus (SLE) patients with deep vein thrombosis (DVT)/ arterial thrombosis (AT) (n = 47), SLE without DVT/AT (n = 60) and healthy donors (HD) (n = 120). OD: optical density.




Figure 3 Comparison of quantitative antibodies to phosphatidylinositol (aPI) IgM (a) and antibodies to aPI IgG (b), antibodies to phosphatidylserine (aPs) IgM (c), aPs IgG (d) analysis detected by line immunoassay (LIA) in systemic lupus erythematosus (SLE) patients with deep vein thrombosis (DVT)/ arterial thrombosis (AT) (n = 47), SLE without DVT/AT (n = 60) and healthy donors (HD) (n = 120). OD: optical density.

Anti-Phospholipid Dot (IgG / IgM)

Why checking for 10 parameters ?

Diagnostic advantages of detection of multiple parameters




Figure 2 Comparison of quantitative anti-cardiolipin (aCL) IgM (a) and anti-beta 2 glycoprotein I (ab2GPI) IgM (b), aCL IgG (c), ab2GPI IgG (d) analysis by line immunoassay (LIA) in systemic lupus erythematosus (SLE) patients with deep vein thrombosis (DVT)/arterial thrombosis (AT) (n=47), SLE without DVT/AT (n=60) and healthy donors (HD) (n=120). OD: optical density.





Figure 3 Comparison of quantitative antibodies to phosphatidylinositol (aPI) IgM (a) and antibodies to aPI IgG (b), antibodies to phosphatidylserine (aPS) IgM (c), aPS IgG (d) analysis detected by line immunoassay (LIA) in systemic lupus erythematosus (SLE) patients with deep vein thrombosis (DVT)/arterial thrombosis (AT) (n=47), SLE without DVT/AT (n=60) and healthy donors (HD) (n=120). OD: optical density.

- Association of non-criteria APL with thrombotic risk in SLE patients
- IgG to PI and PS show significantly higher levels in SLE with deep vein thrombosis/arterial thrombosis than without
- Occurrence of >4 IgG antiphospholipid antibodies is independent risk factor for thrombotic effects

DotDiver Anti-Phospholipid IgG / IgM

DotDiver2.0, Test Strips and Reagents



Antibodies against Phospholipids for the
Diagnosis of the Anti-Phospholipid-Syndrome



Thomas Büttner

Product Management / Customer Support

Medipan GmbH / GA Generic Assays GmbH
Ludwig-Erhard-Ring 3
15827 Blankenfelde-Mahlow OT Dahlewitz
Germany

Phone +49 33708 4417 43
Fax +49 33708 4417 25
E-Mail info@medipan.de



On Track to achieve
IVDR Certification

Anti-Phospholipid-Syndrome

Thank you for your attention

www.medipan.de
www.genericassays.com



Thanks to the company Medi-Lab Kft for
support in the Hungarian Market

Thank you for your attention !