



DiAgam

Specialist in Specific Proteins
Supplier of optimized solution for
liquid clinical chemistry


- **1993 founded by Dr Gilbert Ameryckx as LTD**
- **2003 CE registration for all products**
- **2005 certification for ISO 9001-2000**
- **2005 certification for ISO 13485**
- **2009 certification for ISO 9001-2008**
- **2008 member of the European BIO & Diagnostic Cluster BLOWIN and WalIDX**
- **2009 foundation of nanoQ, LTD sister Company specialized in the synthesis and coating of metallic nanoparticles**
- **2010 implementation of the standard 15185 in the establishment of technical records**

HUMAN RESSOURCES




Qualifications distribution

■ Ph.D.	23%
■ Post Graduates	39%
■ Graduates	38%



Functions distribution

■ Sales and Marketing	34%
■ R+D	22%
■ QA/QC	18%
■ Production	26%



Manufactured Products: specific proteins

Albumin

Alpha1-glycoprotein acid

Alpha1-antitrypsin

Alpha2-macroglobulin

Antithrombin III

Apolipoprotein A1

Apolipoprotein B

ASO

Ceruloplasmin

C3

C4

CRP XL "Wide Range"

CRP XS

Rhumatoid Factor

Ferritin

Fibrinogen

Haptoglobin



IgA, IgG, IgM

Lipoprotein (a) (Lpa)

Microalbumin



Prealbumin

Transferrin



Available for opened analyzers

Advia
Aeroset
Alcyon
Architect
AU series
Cobas Mira
Cobas 6000
Cobas 8000
CX
DXC



Fusion FS 5.1
Hitachi series
Konelab series
Integra series
Modular
Selectra
Targa
Vista
Vitros 5600
...and more...

Meeting LABIMMUN/DIAGAM. 01.12.11.

Strenghts: some exemples

CRP XL
RF
MICROALBUMIN

Strenghts: CRP XL

- **Principle:**

Immunturbidimetry with ready to use reagent, calibrator and control

- **Measuring range:**


0.8 – 300 mg/l

- **Accuracy:**

standardized on ERM-DA474

- **Bar-codes available for:**

Modular, Hitachi, AU, Konelab...



Source: CPTXLHMODE/ACN704/GA/04/04/2008. Gilly

Strenghts: RF

- **Principle:**

Immunoturbidimetry (with colloidal gold particles) with ready to use reagent, calibrator and control

- **Measuring range:**


1.2 – 120 KIU/l

- **Accuracy:**

Standardized on WHO W1066

- **Bar-codes available for:**

Modular, Hitachi, AU, Konelab...



--- : DIAGAM (colloidal gold particles)

--- : Latex based methods

Source: RFCOLHMODE/JOCKER/SQ/03/04/07. BAUDUIN.

Strenghts: MICROALBUMIN

- **Principle:**

Immunturbidimetry with ready to use reagent, calibrator and control

- **Measuring range:**


3.6 – 500 MG/l

- **Accuracy:**

Standardized on ERM-DA470k

- **Bar-codes available for:**

Modular, Hitachi, AU, Konelab...



Source: MATURHMOD ACN005 E/AR/14/04/2009

R+D Products and Technology

CYSTATIN C & D DIMERS

METALLIC
NANOPARTICLES
COATING

