

THE III GENERATION TNI ON PATHFAST® MITSUBISHI : THE SAME QUALITY WITH POCT AS IN THE CENTRAL LABORATORY





F. Manoni, A. Tinello, M. Ercolin, M. Ferrari, E. Giraldo,
V. Temporin, D. Cioetto, C. Vigato, G. Casarotto, P. Hoffer
Laboratory Medicine, Este- Monselice Azienda ULSS 17 Veneto



Summary: Troponin T and I are the reference marker of myocardial injury due to their sensibility and specificity. To make best use of Troponin T and I for clinical purposes the IFCC guidelines recommended a total imprecision equal or less than 10% ($CV \leq 10$) corresponding to a cut-off at 99th percentile of a reference population. These requirements are very critical for many diagnostic system that have high CV at low Troponin concentrations. The aim of this study is to evaluate the PATHFAST® system imprecision around the 99th percentile of the reference population.

Materials and methods : PATHFAST® Mitsubishi performance-distributed by Gepa- is base on CLEIA-MAGTRACTION technology and operates with single use reagent cartridges with an antibody sandwich test format. PATHFAST® Mitsubishi provides 6 results in 17'. According the instructions by the producing company analytical sensibility of the system was evaluated (0.001 ug/L). First of all 123 healthy (no myocardial injury) caucasian subjects (89 male and 34 female) of a medium age of 53 years old (39-66) have been evaluated : corresponding 0.030 $\mu\text{g/L}$ at the 99th percentile. Imprecision was calculated on 11 samples with increasing TnI concentration from 0.001 to 0.066 $\mu\text{g/L}$. The calculation of graphic was executed by Medcalc.




with negative value (0,001-0,020 ug/L)

around cut-off (0,023 - 0,039ug/L)

with feasible positive value (0,042 a 0,066 ug/L)

Result: Test showed an optimal analytical sensibility (0.001 ug/ml) and the imprecision profile is very well covered by the IFCC guidelines at 99th percentile with a $CV \leq 2.1\%$: For all considered TropI concentrations CV was always less than 5.1% ($CV \leq 5.1\%$ see table and graphic).

Samples	Troponin I		Selected samples repeated for 5 times					Mean	DS	CV%
	PATHFAST		run 1	run 2	run 3	run 4	run 5			
1	Negative	0.001	0.001	0.001	0.001	0.001	0.001	0.000	0.00%	
2		0.017	0.016	0.015	0.017	0.016	0.016	0.001	5.16%	
3		0.020	0.020	0.019	0.020	0.021	0.020	0.001	3.54%	
4	Attention Value	0.024	0.023	0.023	0.024	0.022	0.023	0.001	3.61%	
5		0.029	0.030	0.029	0.029	0.030	0.029	0.001	1.86%	
6		0.034	0.034	0.032	0.034	0.033	0.033	0.001	2.68%	
7	Decision Value	0.039	0.039	0.038	0.039	0.040	0.039	0.001	1.81%	
8		0.042	0.043	0.041	0.042	0.040	0.042	0.001	2.74%	
9		0.049	0.051	0.050	0.049	0.050	0.050	0.001	1.68%	
10	Decision Value	0.055	0.056	0.056	0.055	0.057	0.056	0.001	1.50%	
11		0.066	0.063	0.065	0.066	0.067	0.065	0.002	2.32%	



Discussion and conclusions:

From the results it is evident that:

- The excellent performance of analytical sensitivity, the requested precision of TnI at 99° % is beyond the IFCC guideline;
- The excellent analytical sensitivity and precision grant a very low cut-off, essential condition for a better diagnostic sensitivity
- Test is quickly and easy to use with reduced dimension, cartridges "all in one" without any additional preparation; this results makes PATHFAST® particularly suitable as well as for centrallab use as for POCT use.